



EVDS

EQUINE VETERINARY & DENTAL SERVICES PTY LTD

Dr Oliver Liyou & Associates
BVSc (Hons) MACVSc (EqDent)

P | (02) 6642 4700
E | vets@evds.net.au
www.evds.vet

PLEASE SHIP EQUIPMENT TO:
160 Old Lilypool Rd
South Grafton NSW 2460

X RAY GENERATOR REPAIR/SERVICE AUTHORISATION

Please complete and include a copy with your X Ray Generator

Date _____ Vet Name _____

Clinic Name _____ Billing Name (if required) _____

Shipping Address _____

Contact Number _____ Mobile (authorised person) _____

Email _____

X Ray Repair/Service Information

Generator Serial No _____

Model AirRay Mini ☐ AirRay 20 ☐ Vet20-BT ☐

Purchase date (approx) _____

Have you double checked the settings are correct? Yes ☐ No ☐

If Yes, please describe issues _____

Please allow up to 3 months for Services/Repairs + Freight Time

Would you like a Rental Unit sent to you prior to sending your unit? ☐ Yes ☐ No

OFFICE USE ONLY

Date Received at EVDS _____ Parts received _____

Notes _____

Date Returned to Client _____ Freight Tracking _____

AIRRAY MINI GENERATOR

QUICK TROUBLESHOOTING GUIDE

AirRay Mini Default Settings

Set Up → Auto Power-Off 10 min | Screen Lock OFF | Power Save ON (5 min).

Utility → Distance Monitoring OFF | Force Exposure ON | DAP OFF.

SID → cm selected | SID Auto DISABLED | SID Auto-limit OFF.

Lamp → Ready ON for 30 seconds.